




# Alliance 2690/5

Jade C. Byrd  
Waters Corporation



- Solvent Management
  - Quaternary solvent low pressure mixing
  - Vacuum degassing
  - Serial flow path configuration
  - Two inlet check valves with First-in, First-out solvent delivery
  - Programmable piston volume delivering optimum performance without changing pump heads
  - Plunger seal wash
  - Smooth, pulse-free solvent flow




Gradient Proportioning Valve (GPV)


# Low Pressure Mixing

Waters

THE SCIENCE OF WHAT'S POSSIBLE™





**All solvent lines should be primed... even when not in use.**



# Solvent Management: Low Pressure Mixing

Waters


THE SCIENCE OF WHAT'S POSSIBLE™




# Solvent Management vs. Pump

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™

## Solvent Management




## Pump



©2007 Waters Corporation 7

# Serial Pump Design


Waters  
THE SCIENCE OF WHAT'S POSSIBLE™



2690/5 Solvent Management Delivery Theory

## Proof of Performance

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™



5 ng of the five peptide test mixture on column can be easily distinguished from the baseline in the chromatograms generated on a Waters Alliance HPLC.


At 5 ng, the peak is lost in the baseline ripple present in the chromatograms generated on a competitive "traditional design" Binary pump and Quaternary gradient system.

©2007 Waters Corporation 9

## 2695 Sample Management

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™

120 x 2mL vials in 5 x 24 vial carousels




©2007 Waters Corporation 10

## Needle Option

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™



- 2695 utilizes side port needle (B)
  - Stainless steel
  - Standard on 2695
- PerformancePLUS needle
  - Nickel/Cobalt alloy
  - Part Number 700001247



©2007 Waters Corporation 11

## Injector

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™




©2007 Waters Corporation 12

# Syringe Options

Waters

THE SCIENCE OF WHAT'S POSSIBLE™



| Part Number | Syringe (µL) | Steps per µL | Minimum Sample size (µL) | Minimum draw rate µL/min | Maximum Draw rate µL/min |
|-------------|--------------|--------------|--------------------------|--------------------------|--------------------------|
| WAT077343   | 25           | 270          | 0.1 µL                   | 0.1                      | 0.5                      |
| WAT044886   | 250          | 27           | 1.0 µL                   | 1.0                      | 5.0                      |
| WAT077342   | 2500         | 2.8          | 10.0 µL                  | 10.0                     | 50.0                     |

# Automate Priming

Waters

THE SCIENCE OF WHAT'S POSSIBLE™

- Ease of Use - Automation benefits
  - build unattended daily startup routines into sample methods with their Alliance systems – frees up valuable lab time

Sample Set Method: AutoStart gradient 1

| Vial | SampleName       | Inj Vol (µL) | # of Inj | Function         | Method Set / Report Method | Run Time (Minutes) | Next Inj Delay (Minutes) | Auto Additions | analytical_column |
|------|------------------|--------------|----------|------------------|----------------------------|--------------------|--------------------------|----------------|-------------------|
| 1    |                  |              |          | Wet Prime        | Wet Prime A                | 0.20               | 7.50                     |                |                   |
| 2    |                  |              |          | Wet Prime        | Wet Prime B                | 0.20               | 7.50                     |                |                   |
| 3    |                  |              |          | Wet Prime        | Wet Prime C                | 0.20               | 7.50                     |                |                   |
| 4    |                  |              |          | Wet Prime        | Wet Prime D                | 0.20               | 7.50                     |                |                   |
| 5    |                  |              |          | Equilibrate      | equilibrate at 0 flow      | 5.00               | 0.00                     |                |                   |
| 6    |                  |              |          | Equilibrate      | gradient method 1          | 9.00               | 0.00                     |                |                   |
| 7    |                  |              |          | Purge Inj        | gradient method 1          | 7.50               | 0.00                     |                |                   |
| 8    | std check sample | 10.0         | 1        | Inject Standards | gradient method 1          | 10.00              | 6.50                     |                |                   |
| 9    | std check sample | 10.0         | 1        | Inject Standards | gradient method 1          | 10.00              | 6.50                     |                |                   |
| 10   | std check sample | 10.0         | 1        | Inject Standards | gradient method 1          | 10.00              | 6.50                     |                |                   |
| 11   | std check sample | 10.0         | 1        | Inject Standards | gradient method 1          | 10.00              | 6.50                     |                |                   |
| 12   |                  |              |          | Report           | Component Summary          |                    |                          |                |                   |

Direct Functions dialog box:

- 1 Dry Prime
- 2 System Prep
- 3 Condition Column
- 4 Inject Samples
- 5 Wet Prime
- 6 Refresh Syringe
- 7 Wash Needle

## Good Practices

Keeping your LC (and chromatography) working well

## Solvent Preparation

- Mobile phase difficulties are estimated to account for 70% or more of all liquid chromatographic problems
- Use only clean, HPLC grade solvents
- All buffers should be added, pH adjustments completed, then the solvents should be filtered using a 0.45 micron filter
- Use water only from a high quality water purification system, or HPLC grade water
- Do not "top off" bottles, replace bottles when empty
- Pure water tends to grow bacteria which will clog the system and foul the degasser
  - If the system is idle, replace pure water with water that has 10% organic solvent added to inhibit bacterial growth



## Degassing

Waters

THE SCIENCE OF WHAT'S POSSIBLE™

- Regardless of the grade of solvent purchased, bottled reagents contain some particulate matter and dissolved gasses
- Most solvent manufacturers sparge their solvents with nitrogen just prior to bottling
- A freshly opened bottle of solvent is often supersaturated with nitrogen
  - Freshly opened solvents **must** be degassed before use
- Vacuum filtration of all solvents through a .45 micron filter will remove particulates, and degas the solvent

©2007 Waters Corporation 17

## Performing Quality Separations

Waters

THE SCIENCE OF WHAT'S POSSIBLE™



Instrument A --  
**less Band spreading --**  
**(70 $\mu$ L)**  
**MORE PLATES (12,000)**

Instrument B --  
**more Band spreading --**  
**(130 $\mu$ L)**  
**LESS PLATES (7,000)**

©2007 Waters Corporation 18

## What contributes to loss of N?


Waters  
THE SCIENCE OF WHAT'S POSSIBLE™

- Column
  - Void
  - Plugging
- Extra-Column (Non - Column) Band spreading
  - Injection Volume
  - Injector (seal problem)
  - Detector Volume
  - Connecting Tubing
    - Injector to Column
    - Column to Detector
  - End-fittings and Frits

©2007 Waters Corporation 19

## Extra column bandspreading

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™




Tubing was cut poorly at an angle rather than smooth and flush  
Potential for system band spreading due to mixing chamber effect  
- Use tubing cutter

©2007 Waters Corporation 20

**Proper Connection –  
Tubing Flush  
with End-fitting  
and Good Seal**


**Column tubing not  
seated properly --  
creates a void  
(Parker ferrule  
position into a Waters'  
0.130" End-fitting)**



**(PEEK Slip Connectors Easier to Use)**

**Column Installation and Connections**


**Waters Spherisorb®,  
ACQUITY UPLC™  
and many other brands**



**If you are using a stainless steel ferrule on  
stainless tubing, once it is set, it will NOT move**

## Effect of Connecting Tubing Volume (ID and Length) on System Band Spreading

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™



## Tubing Considerations

Waters  
THE SCIENCE OF WHAT'S POSSIBLE™

- Stainless Steel
  - High pressures (w/ SS fittings)
  - Hydrophobic compounds
    - Can "passivate"—acid wash
- PEEK Tubing

Solvents to Avoid with PEEK Tubing

**Attacks PEEK**

Conc.  $\text{HNO}_3$   
Conc.  $\text{H}_2\text{SO}_4$

**Swells PEEK**

DMSO  
THF  
 $\text{CH}_2\text{Cl}_2$

## Reducing Contamination/Ghost Peaks/Carryover

Waters

THE SCIENCE OF WHAT'S POSSIBLE™

- Test pH using an aliquot and then dispose
  - pH probe can contaminate the bottle of mobile phase
- Change your needle wash filter
  - Especially needed if the needle wash has run dry
  - Remember to pre-moisten
- Ensure your needle wash is compatible/good for your sample
  - E.g. salts need some water, maybe acid
- Put a wash step in your gradient
  - Remember the inside of the needle is washed by the gradient
  - Cleans out the column also
- Acid-clean/passivate the HPLC system
  - $\text{H}_3\text{PO}_4$ , Nitric
- Use certified/pre-cleaned vials

©2007 Waters Corporation 25

## Vials

Waters

THE SCIENCE OF WHAT'S POSSIBLE™

- Unreactive Glass
  - Type I borosilicate
    - Least reactive
  - USP
- Consistent dimensions
  - Don't crash needles
- Caps/Septa
  - Core easily
  - Won't detach
  - Release vacuum
    - Pre-slit
  - Low Contaminants
    - PTFE
    - QC tested



26

- Online Library
  - Technical notes
  - Methods and application-specific solutions
- Call/email Waters Home Office
  - Nature of symptom
  - Instrument SN
  - Method parameters
    - Flow rate, solvents/mobile phases, detector wavelength, type of column
  - Operating Pressure
  - Sample type, matrix
- Service Engineer/Plan
  - Plan cost is 85% cost of parts

[www.waters.com/techassist](http://www.waters.com/techassist)

1.800.252.HPLC



Questions?

## What is a Compression Check?

- Automated routine to check for leaks in
  - Sample loop
  - Syringe
  - Needle Seal Pack
  - Valves
- Also detects bubbles in the sample withdrawal path
  - Bubbles will cause problems with system accuracy and precision

## What happens during a Compression Check?

- Sample Loop is isolated from the mobile phase flow
- Sample Loop is vented to atmosphere
- A closed system is created that includes the Sample Loop, Valves, Syringe, and Needle Seal Pack
- The closed system is pressurized to a set pressure by the syringe
  - The volume of solvent required to pressurize the closed system is recorded
- A leak down test is run for a specified time, and the loss of pressure is recorded

Diagnostics  
<Default> Idle

Serial Nu  
Transc  
System  
Primary  
Sample

Compression Check  
Idle

Initializing

Results:

Compressed Volume: 0 µL  
Delta Decay Pressure: 0 psi  
Final Pressure: 0 psi

Start Halt Close

- **Compressible Volume**- the volume of solvent required to pressurize the sample loop to a pre set pressure
  - Typically less than 5 ul using a degassed solvent
  - If the value exceeds 5 ul, possible causes are
    - Mobile phase not degassed
    - Bubble in the syringe
    - Leak in the syringe, sample loop, valves, or needle seal pack
  - Repeating the compression check often reduces the compressible volume



## How do I Interpret the Results?

- **Delta Decay Pressure**-change in sample loop pressure after completing the leak down test
  - Typically less than 30% of the compression check pressure
  - If the value is high, the problem is usually a worn syringe
  - If the value is positive, then there is a leak in one of the valves or in the needle seal pack that is allowing high pressure mobile phase to leak into the sample loop
- **Final Pressure**-sample loop pressure after completing the compression check

## Recommended System Startup Checklist

1. See that power is turned on for the Alliance and detector.
2. Turn on degasser, column and sample temps.
3. Check to see there are sufficient solvents for test.
4. Run a wet prime for each solvent you will use. Prime the seal wash and needle wash.
5. Turn on flow with the needed composition.  
**Allow the delta psi to stabilize.**
6. Purge the injector.
7. Program Empower for method, and prepare samples.  
~ Allow approximately 30 min for equilibration.